


Rabbit Anti-Human Axl Mab For Fc - P30530

Item Number: CM0011696



Introduction

-conjugated rabbit monoclonal antibody targeting human AXL (UniProt P30530) validated for flow cytometry. Ideal for studying AXL-mediated signaling, cancer, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human AXL mAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7
Species	Human
GeneID (Human)	558
GeneID	558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-449 of human AXL. (NP_068713.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag