


Rabbit Anti-Human Axl Monoclonal Antibody For Fc - P30530

Item Number: CM0011338



Introduction

-conjugated rabbit monoclonal antibody targeting human AXL receptor tyrosine kinase (UFO/Tyro7). Designed for flow cytometry, validated in human samples. Explore AXL's roles in cancer, immunity, and viral infections.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human AXL mAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7
Species	Human
GeneID (human)	558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-449 of human AXL (NP_068713.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag