


Rabbit Anti-Human Axl Monoclonal Antibody For Flow Cytometry - P30530

Item Number: CM0008790



Introduction

ABflo 488-conjugated rabbit anti-human AXL monoclonal antibody for flow cytometry. Targets human AXL (P30530), a receptor tyrosine kinase involved in cell survival, proliferation, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human AXL mAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7
Species	Human
Gene ID (Human)	558
Gene ID	558
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 33-449 of human AXL (NP_068713.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag