

Rabbit Anti-Human B7-H2/Cd275 Monoclonal Antibody For Flow Cytometry - O75144

Item Number: CM0004027



Introduction

-conjugated rabbit monoclonal antibody against human B7-H2/CD275 (ICOS ligand), designed for flow cytometry. Ideal for T-cell costimulation and immune checkpoint studies, with human cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human B7-H2/CD275 mAb
Remarks/Alias	B7h; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L
Species	Human
GeneID (Human)	23308
GeneID	23308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human B7-H2/CD275 (NP_056074.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75144
Protein Weight	20kDa/33kDa34kDa
Shipping	Ice bag