

Rabbit Anti-Human Bcl6 Mab - P41182

Item Number: CM0004938



Introduction

Anti-BCL6 monoclonal antibody conjugated to , validated for intracellular flow cytometry (FC) on human samples. Ideal for germinal center B-cell and lymphoma research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human BCL6 mAb
Remarks/Alias	BCL6; BCL5; BCL6A; LAZ3; ZBTB27; ZNF51; B-cell CLL/lymphoma 6
Species	Human
GeneID	604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-350 of human BCL6 (NP_001124317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P41182
Protein Weight	72kDa/78kDa
Shipping	Ice bag