

Rabbit Anti-Human Bcl6 Monoclonal Antibody For Flow Cytometry (Intracellular) - P41182

Item Number: CM0011455



Introduction

-conjugated rabbit monoclonal antibody targeting human BCL6, validated for intracellular flow cytometry. Designed for germinal center B-cell and lymphoma research, with specificity for the P41182 protein.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human BCL6 mAb
Remarks/Alias	BCL5; LAZ3; BCL6A; ZNF51; ZBTB27
Species	Human
GeneID	604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-350 of human/mouse BCL6 (NP_001124317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P41182
Protein Weight	79kDa
Shipping	Ice bag