

Rabbit Anti-Human Bcl6 Monoclonal Antibody For Intracellular Flow Cytometry - P41182

Item Number: CM0011162



Introduction

-conjugated rabbit monoclonal antibody against human BCL6 (SWISS P41182) validated for intracellular flow cytometry. Ideal for studying germinal center biology and transcriptional repression. Supplied by CamelBio for IVD raw material procurement.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human BCL6 mAb
Remarks/Alias	BCL5; LAZ3; BCL6A; ZNF51; ZBTB27
Species	Human
Gene ID (Human)	604
Gene ID	604
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 235-350 of human/mouse BCL6 (NP_001124317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P41182
Protein Weight	79kDa
Shipping	Ice bag