

Rabbit Anti-Human Btn3A1/2/3 Mab - O00481 / P78410 / O00478

Item Number: CM0011732



Introduction

High-quality Rabbit anti-Human BTN3A1/2/3 monoclonal antibody for flow cytometry (FC). Recognizes human BTN3A1/2/3 (CD277), crucial in T-cell activation and immune response. Ideal for immunology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human BTN3A1/2/3 mAb
Remarks/Alias	BTF5; BT3.1; CD277; BTN3.1
Species	Human
GeneID	11119/11118/10384
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 470-630 of human BTN3A1 (NP_008979.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O00481/P78410/O00478
Protein Weight	58kDa
Shipping	Ice bag