

Abflo 488 Rabbit Anti-Human C5Ar1/Cd88 Monoclonal Antibody For Flow Cytometry - P21730

Item Number: CM0004526



Introduction

ABflo 488-conjugated rabbit monoclonal antibody targeting human C5AR1/CD88 (C5aR) validated for flow cytometry. ~39 kDa, human reactivity, shipped on ice. For immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human C5AR1/CD88 mAb
Remarks/Alias	C5A; C5AR; C5R1; CD88
Species	Human
Gene ID (Human)	728
Gene ID	728
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human C5AR1/CD88 (NP_001727.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21730
Protein Weight	39kDa
Shipping	Ice bag