


Rabbit Anti-Human Calreticulin Monoclonal Antibody For Fc (Intra) - P27797

Item Number: CM0011279



Introduction

High-quality -conjugated rabbit monoclonal antibody against human Calreticulin, suitable for intracellular flow cytometry (FC intra). Validated for human samples. Ideal for studying ER chaperone function, protein quality control, and calcium homeostasis. Sourced via CamelBio for IVD research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Calreticulin mAb
Remarks/Alias	RO; CRT; SSA; cC1qR; HEL-S-99n
Species	Human
GeneID (Human)	811
GeneID	811
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P27797
Protein Weight	48kDa
Shipping	Ice bag