


Rabbit Anti-Human Calreticulin Monoclonal Antibody For Flow Cytometry - P27797

Item Number: CM0004016



Introduction

-conjugated rabbit monoclonal antibody targeting human calreticulin for intracellular flow cytometry. Suitable for detecting calreticulin, an ER chaperone involved in protein folding and calcium signaling, in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Calreticulin mAb
Remarks/Alias	RO; CRT; SSA; cC1qR; HEL-S-99n
Species	Human
Gene ID (Human)	811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 335-417 of human Calreticulin (NP_004334.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular)
Cross Reactivity	Human
SWISS	P27797
Protein Weight	48kDa
Shipping	Ice bag