


Rabbit Anti-Human Car/Cxadr Mab For Fc - P78310

Item Number: CM0011508



Introduction

Rabbit anti-Human CAR/CXADR monoclonal antibody for Flow Cytometry (FC). Recognizes human CXADR (CAR), a 40 kDa coxsackievirus and adenovirus receptor involved in cell adhesion and viral entry.

[Learn More](#)

Product Name	Rabbit anti-Human CAR/CXADR mAb
Remarks/Alias	CAR; HCAR; CAR4/6
Species	Human
Gene ID (Human)	1525
Gene ID	1525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-237 of human CAR/CXADR(NP_001329.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78310
Protein Weight	40kDa
Shipping	Ice bag