


Rabbit Anti-Human Car/Cxadr Mab - P78310

Item Number: CM0008315



Introduction

-conjugated rabbit monoclonal antibody targeting human CAR/CXADR (P78310). Validated for flow cytometry. CAR functions as a cell adhesion molecule and viral receptor.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CAR/CXADR mAb
Remarks/Alias	CAR; HCAR; CAR4/6
Species	Human
GeneID (Human)	1525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-237 of human CAR/CXADR (NP_001329.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78310
Protein Weight	40kDa
Shipping	Ice bag