


Rabbit Anti-Human Car/Cxadr Monoclonal Antibody For Flow Cytometry - P78310

Item Number: CM0007957



Introduction

-conjugated rabbit monoclonal antibody targeting human coxsackievirus and adenovirus receptor (CAR/CXADR), validated for flow cytometry. Recognizes an epitope in the extracellular domain (aa20-237). Ideal for research on viral entry, cell adhesion, and epithelial junctions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CAR/CXADR mAb
Remarks/Alias	CAR; HCAR; CAR4/6
Species	Human
GeneID (Human)	1525
GeneID	1525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-237 of human CAR/CXADR(NP_001329.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78310
Protein Weight	40kDa
Shipping	Ice bag