

Rabbit Anti-Human Casr Monoclonal Antibody For Flow Cytometry - P41180

Item Number: CM0009552



Introduction

Rabbit monoclonal -conjugated Anti-Human CaSR antibody for flow cytometry. Validated in human samples; detects the extracellular calcium-sensing receptor, a key regulator of calcium homeostasis and parathyroid function. Ideal for renal and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CaSR mAb
Remarks/Alias	CAR; FHH; FIH; HHC; EIG8; HHC1; NSHPT; PCAR1; hCasR; GPRC2A; HYPOC1
Species	Human
GeneID (Human)	846
GeneID	846
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 211-310 of human CaSR (NP_000379.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P41180
Protein Weight	121kDa
Shipping	Ice bag