

Rabbit Anti-Human Ccd3 Mab For Flow Cytometry (Intracellular) - P07766

Item Number: CM0010383



Introduction

-conjugated rabbit monoclonal antibody against human CD3E (CD3 epsilon) for intracellular flow cytometry. Detects a 23kDa T-cell membrane protein critical in TCR signaling and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human cCD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID	916
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human CD3E (NP_000724.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag