

Rabbit Anti-Human Cd100/Sema4D Mab - Q92854

Item Number: CM0005994



Introduction

Rabbit monoclonal antibody conjugated to , targeting human CD100/SEMA4D. Suitable for immunofluorescence and flow cytometry. Useful for research in semaphorin signaling, immune cell aggregation, and neuronal migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD100/SEMA4D mAb
Remarks/Alias	SEMA4D; C9orf164; CD100; M-sema-G; SEMAJ; coll-4; semaphorin-4D
Species	Human
GeneID (Human)	10507
GeneID	10507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-734 of human Semaphorin-4D/CD100 (NP_006369.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q92854
Protein Weight	82kDa/96kDa
Shipping	Ice bag