


Rabbit Anti-Human Cd100/Sema4D Monoclonal Antibody For If/Icc, Fc - Q92854

Item Number: CM0004358



Introduction

-conjugated rabbit monoclonal antibody against human CD100/SEMA4D, validated for immunofluorescence (IF/ICC) and flow cytometry (FC). Suitable for studying semaphorin-4D receptor signaling, neurite development, and immune cell aggregation. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD100/SEMA4D mAb
Remarks/Alias	SEMA4D; C9orf164; CD100; M-sema-G; SEMAJ; coll-4; semaphorin-4D
Species	Human
Gene ID (Human)	10507
Gene ID	10507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-734 of human Semaphorin-4D/CD100 (NP_006369.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q92854
Protein Weight	82kDa/96kDa
Shipping	Ice bag