

Rabbit Anti-Human Cd101/Igsf2 Monoclonal Antibody For If/Icc And Fc - Q93033

Item Number: CM0004441

Introduction

This rabbit monoclonal antibody targets human CD101/IGSF2 (V7, EWI-101), a T-cell proliferation inhibitor. Suitable for IF/ICC and flow cytometry in immunological studies and assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD101/IGSF2 mAb
Remarks/Alias	V7, IGSF2, EWI-101
Species	Human
Gene ID (Human)	9398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-410 of human CD101/IGSF2 (NP_004249.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q93033
Protein Weight	115kDa
Shipping	Ice bag