

Rabbit Anti-Human Cd101/Igsf2 Monoclonal Antibody For If/Icc, Fc - Q93033

Item Number: CM0011952

Introduction

High-quality -labeled rabbit monoclonal antibody targeting human CD101 (IGSF2, V7) for immunofluorescence and flow cytometry applications. Suitable for T-cell proliferation inhibition research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD101/IGSF2 mAb
Remarks/Alias	V7; IGSF2; EWI-101
Species	Human
GeneID (Human)	9398
GeneID	9398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-410 of human CD101/IGSF2 (NP_004249.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q93033
Protein Weight	115kDa
Shipping	Ice bag