

## Abflo 594 Rabbit Anti-Human Cd104/Integrin Beta-4 Monoclonal Antibody For If/Icc,Fc - P16144

Item Number: CM0012631



### Introduction

High-quality ABflo 594-conjugated rabbit monoclonal antibody against human CD104/Integrin  $\beta$ 4 (ITGB4), validated for immunofluorescence and flow cytometry. Detects a 202kDa protein, cross-reactive with human, ideal for epithelial cell and hemidesmosome research.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD104/Integrin $\beta$ 4 mAb                                                                                          |
| Remarks/Alias    | CD104; GP150; JEB5A; JEB5B                                                                                                              |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 3691                                                                                                                                    |
| GeneID           | 3691                                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 611-710AA of human CD104/Integrin $\beta$ 4 (NP_000204.3) |
| Source           | Rabbit                                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                                   |
| Application      | IF/ICC,FC                                                                                                                               |
| Cross Reactivity | Human                                                                                                                                   |
| SWISS            | P16144                                                                                                                                  |
| Protein Weight   | 202kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                 |