

Rabbit Anti-Human Cd104/Integrin B4 Monoclonal Antibody For If/Icc, Fc - P16144

Item Number: CM0010117



Introduction

-conjugated rabbit monoclonal antibody against human Integrin $\beta 4$ (CD104, GP150). Validated for IF/ICC and FC applications. Functions as laminin receptor and component of hemidesmosomes.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD104/Integrin $\beta 4$ mAb
Remarks/Alias	CD104; GP150; JEB5A; JEB5B
Species	Human
GeneID (human)	3691
GeneID	3691
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 611-710AA of human CD104/Integrin $\beta 4$ (NP_000204.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P16144
Protein Weight	202kDa
Shipping	Ice bag