

Rabbit Anti-Human Cd104/Integrin B4 Monoclonal Antibody For If/Icc, Fc - P16144

Item Number: CM0004033



Introduction

This -conjugated rabbit monoclonal antibody targets human CD104/Integrin β 4 (SWISS: P16144), validated for IF/ICC and FC applications. It recognizes the extracellular domain (aa 611-710) of ITGB4, suitable for epithelial cell adhesion and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD104/Integrin β 4 mAb
Remarks/Alias	CD104; GP150; JEB5A; JEB5B
Species	Human
GeneID (Human)	3691
GeneID	3691
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 611-710AA of human CD104/Integrin β 4 (NP_000204.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC,FC
Cross Reactivity	Human
SWISS	P16144
Protein Weight	202kDa
Shipping	Ice bag