

Rabbit Anti-Human Cd105/Endoglin Monoclonal Antibody For Flow Cytometry - P17813

Item Number: CM0004053



Introduction

-conjugated rabbit monoclonal antibody against human CD105/Endoglin, validated for flow cytometry (FC). Endoglin is a transmembrane glycoprotein involved in angiogenesis and TGF-beta signaling. Suitable for endothelial cell and immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD105/Endoglin mAb
Remarks/Alias	ENG; END; HHT1; ORW1; endoglin
Species	human
GeneID	2022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-337 of human CD105/Endoglin (NP_001108225.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P17813
Protein Weight	67kDa/70kDa
Shipping	Ice bag