

Rabbit Anti-Human Cd112/Nectin-2 Mab - Q92692

Item Number: CM0012031



Introduction

-conjugated rabbit monoclonal antibody targeting human CD112/Nectin-2 for flow cytometry. Validated for human samples, this antibody enables T-cell signaling, immune checkpoint, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD112/Nectin-2 mAb
Remarks/Alias	NECTIN2; CD112; HVEB; PRR2; PVRL2; PVRR2; nectin-2
Species	Human
GeneID (Human)	5819
GeneID	5819
Immunogen	Recombinant protein corresponding to a sequence within amino acids 32-158 of human CD112/Nectin-2 (NP_001036189.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q92692
Protein Weight	51kDa/57kDa
Shipping	Ice bag