

Rabbit Anti-Human Cd112/Nectin-2 Monoclonal Antibody For Flow Cytometry - Q92692

Item Number: CM0012243



Introduction

-conjugated rabbit monoclonal antibody targeting human CD112/Nectin-2, validated for flow cytometry. Recognizes an extracellular epitope within amino acids 32-158. Suitable for T-cell signaling and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD112/Nectin-2 mAb
Remarks/Alias	NECTIN2; CD112; HVEB; PRR2; PVRL2; PVRR2; nectin-2
Species	Human
Gene ID (Human)	5819
Gene ID	5819
Immunogen	Recombinant protein corresponding to a sequence within amino acids 32-158 of human CD112/Nectin-2 (NP_001036189.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q92692
Protein Weight	51kDa/57kDa
Shipping	Ice bag