

Rabbit Anti-Human Cd116 Monoclonal Antibody For Fc - P15509

Item Number: CM0011788



Introduction

-conjugated rabbit monoclonal antibody against human CD116 (GM-CSF receptor α). Validated for flow cytometry. Ideal for hematopoietic cell studies and receptor expression analysis. SWISS: P15509.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD116 mAb
Remarks/Alias	GMR; CD116; CSF2R; SMDP4; CDw116; CSF2RX; CSF2RY; GMCSFR; CSF2RAX; CSF2RAY; alphaGMR; GMR-alpha; GMCSFR-alpha; GM-CSF-R-alpha
Species	Human
Gene ID	1438
Immunogen	Recombinant protein of human CD116 (NP_006131.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15509
Protein Weight	46kDa
Shipping	Ice bag