

Rabbit Anti-Human Cd117/C-Kit Mab For Flow Cytometry - P10721

Item Number: CM0004440



Introduction

Rabbit anti-Human CD117/c-Kit mAb validated for flow cytometry (FC). Recognizes human CD117/c-Kit (SWISS P10721), a receptor tyrosine kinase crucial for hematopoiesis and stem cell maintenance. Ideal for hematopoietic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD117/c-Kit mAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC
Species	Human
GeneID (Human)	3815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-520 of human CD117/c-Kit (NP_000213.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag