

Rabbit Anti-Human Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P10721

Item Number: CM0012326



Introduction

Rabbit monoclonal antibody targeting human CD117/c-Kit, conjugated with for flow cytometry. Recognizes recombinant extracellular domain (aa 26-520). Reacts with human cells, ideal for hematopoietic stem cell and mast cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD117/c-Kit mAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC
Species	Human
Gene ID	3815
Immunogen	Recombinant fusion protein containing amino acids 26-520 of human CD117/c-Kit (NP_000213.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag