

Rabbit Anti-Human Cd119/Ifngr1 Monoclonal Antibody For Fc - P15260

Item Number: CM0008019



Introduction

Rabbit anti-Human CD119/IFNGR1 monoclonal antibody for flow cytometry (FC). Specific for human interferon gamma receptor 1, a key immune regulator. Suitable for antiviral, antitumor, and infectious disease studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD119/IFNGR1 mAb
Remarks/Alias	IFNGR1; CD119; IFNGR; IMD27A; IMD27B; interferon gamma receptor 1
Species	Human
Gene ID (Human)	3459
Gene ID	3459
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-245 of human CD119/IFNGR1. (NP_000407.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15260
Protein Weight	21kDa/54kDa
Shipping	Ice bag