

Rabbit Anti-Human Cd119/Ifngr1 Monoclonal Antibody For Flow Cytometry (Fc) - P15260

Item Number: CM0010405



Introduction

-conjugated rabbit monoclonal antibody against human CD119/IFNGR1. Designed for flow cytometry (FC) and validated using recombinant protein immunogen. Research-use only IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD119/IFNGR1 mAb
Remarks/Alias	IFNGR1; CD119; IFNGR; IMD27A; IMD27B; interferon gamma receptor 1
Species	Human
Gene ID	3459
Immunogen	Recombinant protein of human CD119/IFNGR1
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15260
Protein Weight	21kDa/54kDa
Shipping	Ice bag