

Rabbit Anti-Human Cd119/Ifngr1 Monoclonal Antibody For Flow Cytometry - P15260

Item Number: CM0012193



Introduction

-conjugated rabbit anti-human CD119/IFNGR1 monoclonal antibody for flow cytometry (FC). Engineered against extracellular domain AA 18-245; ideal for studying interferon gamma receptor-mediated immune responses. - P15260

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD119/IFNGR1 mAb
Remarks/Alias	CD119; IFNGR; IMD27A; IMD27B
Species	Human
GeneID (Human)	3459
GeneID	3459
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-245 of human CD119/IFNGR1. (NP_000407.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15260
Protein Weight	21kDa/54kDa
Shipping	Ice bag