


Rabbit Anti-Human Cd11A Mab For Fc - P20701

Item Number: CM0004945



Introduction

High-quality -conjugated rabbit anti-human CD11a monoclonal antibody for flow cytometry (FC). Recognizes the LFA-1 alpha chain (integrin alpha-L), a key leukocyte adhesion receptor. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD11a mAb
Remarks/Alias	CD11A; LFA-1; LFA1A
Species	Human
GeneID (Human)	3683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-340 of human CD11a (NP_002200.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20701
Protein Weight	118kDa/129kDa/134kDa
Shipping	Ice bag