


Rabbit Anti-Human Cd11A Mab - P20701

Item Number: CM0004483



Introduction

Rabbit anti-Human CD11a monoclonal antibody for flow cytometry. Targets integrin alpha-L (LFA-1), a leukocyte adhesion receptor involved in immune cell interactions. Reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD11a mAb
Remarks/Alias	CD11A; LFA-1; LFA1A
Species	Human
GeneID (Human)	3683
GeneID	3683
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 26-340 of human CD11a (NP_002200.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20701
Protein Weight	118kDa/129kDa/134kDa
Shipping	Ice bag