


Rabbit Anti-Human Cd120B/Tnfrsf1B Mab For Fc - P20333

Item Number: CM0006095



Introduction

-conjugated rabbit monoclonal antibody targeting human CD120b/TNFRSF1B (TNF-R2/p75). Validated for flow cytometry, suitable for studying TNF signaling, apoptosis, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD120b/TNFRSF1B mAb
Remarks/Alias	p75; TBPII; TNFBR; TNFR2; CD120b; TNFR1B; TNFR80; TNF-R75; p75TNFR; TNF-R-II
Species	Human
GeneID (human)	7133
GeneID	7133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-257 of human CD120b/TNFRSF1B (NP_001057.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20333
Protein Weight	28kDa/48kDa
Shipping	Ice bag