

Rabbit Anti-Human Cd131 Mab For Fc - P32927

Item Number: CM0011672



Introduction

-conjugated rabbit monoclonal antibody targeting human CD131 (CSF2RB/IL3RB/IL5RB), the common beta subunit for IL-3, IL-5, and GM-CSF receptors. Validated for flow cytometry (FC). Ideal for hematopoietic and cytokine signaling research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD131 mAb
Remarks/Alias	CD131; IL3RB; IL5RB; SMDP5; CDw131; betaGMR
Species	Human
GeneID (Human)	1439
GeneID	1439
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-443 of human CD131 (NP_000386.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32927
Protein Weight	97kDa
Shipping	Ice bag