


Rabbit Anti-Human Cd134/Ox40 Mab For Fc - P43489

Item Number: CM0007860



Introduction

-conjugated rabbit monoclonal antibody targeting human CD134/OX40 (TNFRSF4). Designed for flow cytometry (FC) in T-cell activation and immunotherapy research. Specific to human CD134.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD134/OX40 mAb
Remarks/Alias	OX40; ACT35; CD134; IMD16; TXGP1L
Species	Human
GeneID (Human)	7293
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-216 of human CD134/OX40 (NP_003318.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43489
Protein Weight	29kDa
Shipping	Ice bag