

Rabbit Anti-Human Cd138 (Syndecan-1) Monoclonal Antibody For Flow Cytometry - P18827

Item Number: CM0009030



Introduction

Rabbit monoclonal antibody targeting human CD138/Syndecan-1 (SDC1), validated for flow cytometry. Recombinant immunogen (aa 23-251) ensures high specificity for exosome and cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD138 mAb
Remarks/Alias	SDC; CD138; SYND1; syndecan
Species	Human
GeneID (Human)	6382
GeneID	6382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-251 of human CD138 (NP_002988.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag