

Rabbit Anti-Human Cd138/Syndecan-1 Monoclonal Antibody For Fc - P18827

Item Number: CM0011364



Introduction

-conjugated rabbit monoclonal antibody targeting human CD138 (Syndecan-1) for flow cytometry. Specific for human syndecan-1, a transmembrane proteoglycan involved in exosome biogenesis and dental development. Ideal for cancer and exosome research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD138/Syndecan-1 mAb
Remarks/Alias	SDC; CD138; SYND1; syndecan
Species	Human
GeneID (Human)	6382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-251 of human CD138/Syndecan-1 (NP_002988.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag