

## Rabbit Anti-Human Cd13/Anpep Monoclonal Antibody For Flow Cytometry - P15144

Item Number: CM0005667



### Introduction

-conjugated rabbit monoclonal anti-human CD13/ANPEP antibody validated for flow cytometry. CD13/ANPEP is an aminopeptidase involved in peptide digestion and a receptor for HCoV-229E. Ideal for immunology and infectious disease research.

[Learn More](#)

| Parameter        | Value                                                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD13/ANPEP mAb                                                                                                           |
| Remarks/Alias    | APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150                                                                                       |
| Species          | Human                                                                                                                                      |
| GeneID (human)   | 290                                                                                                                                        |
| Immunogen        | Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP (NP_001141.2) |
| Source           | Rabbit                                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                                      |
| Application      | FC                                                                                                                                         |
| Cross Reactivity | Human                                                                                                                                      |
| SWISS            | P15144                                                                                                                                     |
| Protein Weight   | 110kDa                                                                                                                                     |
| Shipping         | Ice bag                                                                                                                                    |