

Rabbit Anti-Human Cd13/Anpep Monoclonal Antibody For Flow Cytometry - P15144

Item Number: CM0005755



Introduction

-conjugated rabbit monoclonal antibody targeting human CD13/ANPEP for flow cytometry. Detects cell surface aminopeptidase N on leukocytes and tumor cells. Hu-reactivity, 110kDa. Reliable for immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD13/ANPEP mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150
Species	Human
GeneID (human)	290
GeneID	290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP (NP_001141.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag