

Rabbit Anti-Human Cd13/Anpep Monoclonal Antibody For Flow Cytometry - P15144

Item Number: CM0009788



Introduction

High-quality -conjugated rabbit monoclonal antibody targeting human CD13/ANPEP, validated for flow cytometry (FC). Recombinant protein immunogen covering amino acids 69-967. Recognizes native conformation on the cell surface. Suitable for immunophenotyping and diagnostics research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD13/ANPEP mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150
Species	Human
GeneID (Human)	290
GeneID	290
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP(NP_001141.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag