

Rabbit Anti-Human Cd13/Anpep Monoclonal Antibody For Flow Cytometry - P15144

Item Number: CM0004909



Introduction

Rabbit anti-Human CD13/ANPEP monoclonal antibody for flow cytometry, targeting the 110 kDa aminopeptidase N (CD13) expressed on human cells. Suitable for research in angiogenesis, cancer, and viral infection.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD13/ANPEP mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150
Species	Human
Gene ID (Human)	290
Gene ID	290
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag