


Rabbit Anti-Human Cd140B/Pdgfr β Monoclonal Antibody For Flow Cytometry - P09619

Item Number: CM0009433



Introduction

Monoclonal antibody conjugated to , validated for flow cytometry detection of human CD140b/PDGFR β . Suitable for studying platelet-derived growth factor receptor beta in cell proliferation, vascular development, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD140b/PDGFR β mAb
Remarks/Alias	PDGFRB; CD140B; IBGC4; IMF1; JTK12; KOGS; PDGFR; PDGFR-1; PDGFR1; PENTT; platelet-derived growth factor receptor beta
Species	Human
GeneID	5159
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 33-530 of human CD140b/PDGFRB(NP_002600.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09619
Protein Weight	37kDa/123kDa
Shipping	Ice bag