


Rabbit Anti-Human Cd141/Thrombomodulin Monoclonal Antibody For Flow Cytometry - P07204

Item Number: CM0008125



Introduction

-conjugated rabbit monoclonal antibody targeting human CD141/Thrombomodulin. Validated for flow cytometry, suitable for hemostasis, coagulation, and inflammation studies on endothelial cells. Research-use diagnostic raw material.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD141/Thrombomodulin mAb
Remarks/Alias	TM; THRM; AHUS6; BDCA3; CD141; BDCA-3; THPH12
Species	Human
Gene ID (Human)	7056
Gene ID	7056
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-515 of human CD141/Thrombomodulin (NP_000352.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07204
Protein Weight	60kDa
Shipping	Ice bag