

Rabbit Anti-Human Cd142 Mab For Flow Cytometry - P13726

Item Number: CM0012492



Introduction

-conjugated rabbit monoclonal antibody targeting human CD142/Tissue Factor, validated for flow cytometry (FC). Detects a 33kDa membrane protein essential for coagulation. Ideal for hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD142 mAb
Remarks/Alias	TF; TFA; CD142
Species	Human
Gene ID (Human)	2152
Gene ID	2152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-251 of human CD142 (NP_001984.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13726
Protein Weight	33kDa
Shipping	Ice bag