


Rabbit Anti-Human Cd142 Monoclonal Antibody For Flow Cytometry - P13726

Item Number: CM0009658



Introduction

-conjugated Rabbit anti-Human CD142 monoclonal antibody for flow cytometry (FC) application. Ideal for coagulation factor III/Tissue factor research; conjugate for direct detection. Human reactivity.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD142 mAb
Remarks/Alias	TF; TFA; CD142
Species	Human
GeneID (Human)	2152
GeneID	2152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-251 of human CD142 (NP_001984.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13726
Protein Weight	33kDa
Shipping	Ice bag