


Rabbit Anti-Human Cd142 (Tissue Factor/Tf) Monoclonal Antibody For Flow Cytometry - P13726

Item Number: CM0008704



Introduction

Rabbit monoclonal antibody against human CD142 (Tissue Factor, TF), conjugated to , validated for Flow Cytometry. Ideal for investigating coagulation, hemostasis, and TF expression in human cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD142 mAb
Remarks/Alias	TF; TFA; CD142
Species	Human
Gene ID (Human)	2152
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 34-251 of human CD142 (NP_001984.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13726
Protein Weight	33kDa
Shipping	Ice bag