

Rabbit Anti-Human Cd144/VE Cadherin Mab For Fc - P33151

Item Number: CM0006046

Introduction

Rabbit monoclonal anti-human VE-cadherin (CD144) antibody, clone 7B4, validated for flow cytometry. Suitable for endothelial cell adhesion studies, vascular biology, and cell junction research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD144/VE Cadherin mAb
Remarks/Alias	7B4; CD144
Species	Human
GenelD (Human)	1003
GenelD	1003
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-593 of human CD144/VE Cadherin (NP_001786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33151
Protein Weight	88kDa
Shipping	Ice bag