

Rabbit Anti-Human Cd144/VE Cadherin Monoclonal Antibody For Fc - P33151

Item Number: CM0007756



Introduction

-conjugated rabbit monoclonal antibody against human CD144/VE-cadherin for flow cytometry. Validated in human samples; suitable for studying endothelial adhesion and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD144/VE Cadherin mAb
Remarks/Alias	7B4; CD144
Species	Human
GeneID (Human)	1003
GeneID	1003
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 48-593 of human CD144/VE Cadherin (NP_001786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33151
Protein Weight	88kDa
Shipping	Ice bag